

Venous Return Cannulae, Single Stage, Right Angle Lighthouse Tip, Wire-reinforced Tubing, 26 Fr, RV-41026



The Venous Return Cannula is indicated for use in drainage of the superior and inferior vena cava during cardiopulmonary bypass surgery. The Catheter consists of a flexible polyvinylchloride tubing wire reinforced tapered to a flow tip. The Catheters are available with a 90-degree bended distal section with a fluted distal flow tip.

| | |
|-----------------------------|----------|
| ORDER CODE | RV-41026 |
| SIZE | 26 Fr |
| AVAILABLE COATED | no |
| QTY PER BOX | 10 |
| LATEX-FREE | yes |
| CONNECTOR ACCEPTANCE | 3/8 |