

# Venous Return Cannulae, Single Stage, Straight Lighthouse Tip, Wire-reinforced Tubing, 12 Fr, RV-40012



The Venous Return Cannula is indicated for use in drainage of the superior and inferior vena cava during cardiopulmonary bypass surgery. The Catheter consists of a flexible polyvinylchloride tubing wire reinforced tapered to a flow tip. The Catheters are available with a 90-degree bended distal section with a fluted distal flow tip.

|                              |          |
|------------------------------|----------|
| <b>ORDER CODE</b>            | RV-40012 |
| <b>SIZE</b>                  | 12 Fr    |
| <b>AVAILABLE COATED</b>      | no       |
| <b>EFFECTIVE LENGTH (CM)</b> | 30,5     |
| <b>QTY PER BOX</b>           | 10       |
| <b>LATEX-FREE</b>            | yes      |
| <b>CONNECTOR ACCEPTANCE</b>  | 1/4      |