

# Venous Return Cannulae, Single Stage, Straight Lighthouse Tip, Wire-reinforced Tubing, 18 Fr, RV-40018



The Venous Return Cannula is indicated for use in drainage of the superior and inferior vena cava during cardiopulmonary bypass surgery. The Catheter consists of a flexible polyvinylchloride tubing wire reinforced tapered to a flow tip. The Catheters are available with a 90-degree bended distal section with a fluted distal flow tip.

Order Code RV-40018

Size 18 Fr

Available Coated no

Effective Length (cm) 36

Qty Per Box 10

Latex-Free yes

Connector Acceptance 1/4"-3/8