

Venous Return Cannulae, Single Stage, Straight Lighthouse Tip, Wire-reinforced Tubing, 20 Fr, RV-40020



The Venous Return Cannula is indicated for use in drainage of the superior and inferior vena cava during cardiopulmonary bypass surgery. The Catheter consists of a flexible polyvinylchloride tubing wire reinforced tapered to a flow tip. The Catheters are available with a 90-degree bended distal section with a fluted distal flow tip.

| | |
|------------------------------|----------|
| ORDER CODE | RV-40020 |
| SIZE | 20 Fr |
| AVAILABLE COATED | no |
| EFFECTIVE LENGTH (CM) | 36 |
| QTY PER BOX | 10 |
| LATEX-FREE | yes |
| CONNECTOR ACCEPTANCE | 1/4"-3/8 |