

# Venous Return Cannulae, Single Stage, Straight Lighthouse Tip, Wire-reinforced Tubing, 32 Fr, RV-40032



The Venous Return Cannula is indicated for use in drainage of the superior and inferior vena cava during cardiopulmonary bypass surgery. The Catheter consists of a flexible polyvinylchloride tubing wire reinforced tapered to a flow tip. The Catheters are available with a 90-degree bended distal section with a fluted distal flow tip.

Order Code                      RV-40032

Size                                32 Fr

Available Coated              no

Effective Length (cm)      41

Qty Per Box                    10

Latex-Free                      yes

Connector Acceptance      3/8 in