

Vents, PVC Open Blunt Tip, with removeable Stylet, with Two Positions Latching Device, 12 Fr, VT-31212



The Vent Catheters are indicated for use during cardiopulmonary bypass surgery to vent/drain the left ventricle.

|                       |          |
|-----------------------|----------|
| ORDER CODE            | VT-31212 |
| SIZE                  | 12 Fr    |
| AVAILABLE COATED      | no       |
| EFFECTIVE LENGTH (CM) | 27,2     |
| QTY PER BOX           | 10       |
| LATEX-FREE            | yes      |