

Vents, PVC Open Blunt Tip, with removeable Stylet, with Two Positions Latching Device, 18 Fr, VT-53218



The Vent Catheters are indicated for use during cardiopulmonary bypass surgery to vent/drain the left ventricle.

ORDER CODE	VT-53218
SIZE	18 Fr
AVAILABLE COATED	no
EFFECTIVE LENGTH (CM)	36
QTY PER BOX	10
LATEX-FREE	yes
NO DOP	yes



Contact us Complete the form to get in touch with a LivaNova representative.

Contact Us