


# Vents, PVC Right Angle, Bullet Tip, without Stylet, 10 Fr, VT-84510



The Vent Catheters are indicated for use during cardiopulmonary bypass surgery to vent/drain the left ventricle.

ORDER CODE	VT-84510
SIZE	10 Fr
AVAILABLE COATED	no
EFFECTIVE LENGTH (CM)	31,8
QTY PER BOX	10
LATEX-FREE	yes
NO DOP	yes



**Contact us** Complete the form to get in touch with a LivaNova representative.

Contact Us