


Vents, PVC Straight, Bullet Tip, with Guidewire Stylet, 13 Fr, VT-84413



The Vent Catheters are indicated for use during cardiopulmonary bypass surgery to vent/drain the left ventricle.

ORDER CODE	VT-84413
SIZE	13 Fr
AVAILABLE COATED	no
EFFECTIVE LENGTH (CM)	31,8
QTY PER BOX	10
LATEX-FREE	yes



Contact us Complete the form to get in touch with a LivaNova representative.

Contact Us