

Arterial Cannulae, 135 Curved Tip, with Suture Collar, Wire-reinforced Tubing, 5 Fr/1.5 mm, A272-15N



The Pediatric cannulae are intended to be used in both the venous and the arterial line of the extracorporeal circuit during cardiopulmonary surgery for periods of up to six hours. Specifically, Arterial Cannulae are designed to return blood to the patient during cardiopulmonary bypass procedures. Depending on country registration, an optional phospholipid surface treatment may be available. Note on ordering: coated products have a “B” before their reference code. Please, contact your local sales representative.

ORDER CODE	A272-15N
SIZE	5 Fr/1.5 mm
AVAILABLE COATED	Depending on country registration
EFFECTIVE LENGTH (CM)	22
QTY PER BOX	10
LATEX-FREE	yes



Contact us Complete the form to get in touch with a LivaNova representative.

Contact Us