

Arterial Cannulae, Curved Tip, with Suture Collar, Wire-reinforced Tubing, 9 Fr/3.0 mm, A222-30B



The Pediatric cannulae are intended to be used in both the venous and the arterial line of the extracorporeal circuit during cardiopulmonary surgery for periods of up to six hours. Specifically, Arterial Cannulae are designed to return blood to the patient during cardiopulmonary bypass procedures. Depending on country registration, an optional phospholipid surface treatment may be available. Note on ordering: coated products have a “B” before their reference code. Please, contact your local sales representative.

ORDER CODE	A222-30B
SIZE	9 Fr/3.0 mm
AVAILABLE COATED	Depending on country registration
EFFECTIVE LENGTH (CM)	22
QTY PER BOX	10
LATEX-FREE	yes



Contact us Complete the form to get in touch with a LivaNova representative.

Contact Us