

Suction Products, Fluted Tip, Plastic Handle, 8 Fr, S354-25



Livanova Pediatric Suckers are intended for intra-cardiac pediatric usage to ensure better visualization of the surgical field and are available with ergonomic plastic grip or rigid tubing. Depending on country registration, an optional phospholipid coating may be available to help improve biocompatibility. Note on ordering: coated products have a “B” before their reference code. Please, contact your local sales representative.

ORDER CODE	S354-25
SIZE	8 Fr
AVAILABLE COATED	Depending on country registration
QTY PER BOX	10
LATEX-FREE	yes
CONNECTOR	1/4



Contact us Complete the form to get in touch with a LivaNova representative.

Contact Us