

Suction Products, Suction Wands, Rigid Fluted Tip with Metal Tube and Plastic Handle, 12.8 Fr, SU-10101



Livanova Rigid- and Flex-Tip Suction Wands are intended for intra-cardiac use. They are made of metal suction tubes with plastic handles. The Low Profile Suction Wand has a malleable low-profile metal suction tube with a plastic handle, pre-attached tubing and clamp, and terminates with a slip-on 1/4" connector. Depending on country registration, an optional phospholipid coating may be available to help improve biocompatibility. Note on ordering: coated products have a “B” before their reference code. Please, contact your local sales representative.

ORDER CODE	SU-10101
SIZE	12.8 Fr
AVAILABLE COATED	no
QTY PER BOX	10
LATEX-FREE	yes



Contact us Complete the form to get in touch with a LivaNova representative.

Contact Us