

Venous Return Cannulae, Single Stage, Bullet Tip Venous Return Cannulae, Wire-reinforced Tubing, 18 Fr, V122- 18



The Pediatric cannulae are intended to be used in both the venous and the arterial line of the extracorporeal circuit during cardiopulmonary surgery for periods of up to six hours. Specifically, Venous Cannulae are designed to drain blood from the patient during cardiopulmonary bypass procedures. Depending on country registration, an optional phospholipid surface treatment may be available. Note on ordering: coated products have a “B” before their reference code. Please, contact your local sales representative.

Order Code V122-18

Qty Per Box 10

Latex-Free yes