


# Venous Return Cannulae, Single Stage, Right Angle Lighthouse Tip, Wire-reinforced Tubing, 16 Fr, RV-41016



The Venous Return Cannula is indicated for use in drainage of the superior and inferior vena cava during cardiopulmonary bypass surgery. The Catheter consists of a flexible polyvinylchloride tubing wire reinforced tapered to a flow tip. The Catheters are available with a 90-degree bended distal section with a fluted distal flow tip.

ORDER CODE	RV-41016
SIZE	16 Fr
AVAILABLE COATED	no
QTY PER BOX	10
LATEX-FREE	yes
CONNECTOR ACCEPTANCE	1/4



**Contact us** Complete the form to get in touch with a LivaNova representative.

Contact Us